

**CERTIFICATE OF ANALYSIS**

## Canine Distemper Virus (CDV)

FA Control Slide

Catalog No.:	SLD-FAC-CDV
Lot:	P171006-001
Expiration:	21 May 2025

**Quality Control Method:** Direct FA on CDV polyclonal conjugate (catalog no. CJ-F-CDV-1ML or 10ML) and CDV monoclonal conjugate (catalog no. CJ-F-CDV-MAB-1ML or 10ML).

**Specific Reaction:** 2–4+ signal on the positive well and negative on the negative well with no background. There are 0 - 50 infected cells per high-power field.

**Other Comments:** NA

**QC Manager:**



David Hyatt

**Pattern Of Fluorescence:** Individual cells with inclusion bodies and some plasma membrane fluorescence.

Technical Data Sheet Version: 2

**FOR IN VITRO LABORATORY USE ONLY.**

H:\Master Documents\COA (VMRD)\2-well Slides\CDV\SLD-FAC-CDV\_P171006-001\_250521.docx

01 June 2021