

CERTIFICATE OF ANALYSIS

# Canine Parvovirus (CPV)


FA Control Slide

|              |             |
|--------------|-------------|
| Catalog No.: | SLD-FAC-CPV |
| Lot:         | P160524-001 |
| Expiration:  | 22 May 2024 |

Quality Control Method: Direct FA using CPV direct conjugate (catalog no. CJ-F-CPV-MAB-10ML).

Specific Reaction: 2-4+ signal on the positive well, negative on the negative well with trace background. There are greater than 25 infected cells per high powered field.

Other Comments: NA

QC Manager:   
David Hyatt

Pattern of Fluorescence: Single cells with cytoplasmic and nuclear fluorescence.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\2-well Slides\CPV\SLD-FAC-CPV\_P160524-001\_240522.docx

28 May 2020