

CERTIFICATE OF ANALYSIS

Rabies Recombinant Nucleoprotein (rNP)

FA Control Slide

Catalog No.:	SLD-FAC-RAB
Lot:	P170424-003
Expiration:	18 May 2021

Quality Control Method: Direct FA using Rabies FITC in-house conjugate.

Specific Reaction: 2-4+ signal on the positive well, negative on the negative well with no background. There are 0-20 infected cells per high powered field.

Other Comments: NA

QC Manager: *Will Thurman* Will Thurman

QC Release Date: 19 May 2021

Pattern of Fluorescence: Smooth and particulate cytoplasmic fluorescence of single and multinucleated cells.

Technical Data Sheet Version: Version #3

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\2-well Slides\Rabies\SLD-FAC-RAB_P170424-003_210518.docx

26 June 2017