

CERTIFICATE OF ANALYSIS

Canine Coronavirus (CCV)

FA Substrate Slide

Catalog No.:	SLD-IFA-CCV
Lot:	P150402-009
Expiration:	14 May 2023

Quality Control Method: Indirect FA using CCV antiserum (catalog no PAB-CCV) and anti-porcine affinity purified IgG FITC conjugate (catalog no CJ-F-PORG-AP-1ML or 10ML).

Specific Reaction: 1 - 3+ signal at 1/400, no background. There are 0 to 50 infected cells or syncytia per high-power field.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 14 May 2019

Pattern Of Fluorescence: Large fluorescent syncytia with some single cells in between.

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.

16 May 2019