

**VMRD**

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: vmrd@vmrd.comWeb site: <http://www.vmrd.com>

Certificate of Analysis

Canine Parainfluenza Virus (CPI)

Anti-Viral Antiserum

CATALOG NO.: 210-70-CPI

VOLUME: 2 ml

LOT: P080709-005

EXPIRATION DATE: 24 July 2010

VIRUS: Canine Parainfluenza Virus (CPI)

DESCRIPTION: Anti-CPI polyclonal antiserum. Liquid. Porcine origin.


INTENDED USE: Useful for IFA. Not suitable for cell culture serum neutralization because it contains azide.

QUALITY CONTROL METHOD: Indirect FA using VMRD, Inc. CPI FA 10-well Slides (catalog no. SLD-IFA-CPI); Anti-Pig IgG FITC Conjugate (catalog no. 085-10).

Specific Reaction: 4+ positive fluorescence at 1:3200. Endpoint at 1:12,800.

Other Reactions or Comments: <1:50 on CPV, CDV, CHV. Endpoint @ 1:200 on CAV.

STORAGE: This antiserum is provided in liquid form and should be stored at 4-8°C. DO NOT FREEZE! If antiserum becomes cloudy, it should be discarded. This antiserum contains 0.09% sodium azide as a preservative.



FOR *IN VITRO* LABORATORY USE ONLY.

WARRANTY: VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.