



VMRD

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: vmrd@vmrd.com

Web site: <http://www.vmrd.com>

Certificate of Analysis

Porcine Circovirus (PCV2)

Anti-Viral Antiserum


CATALOG NO.: 210-70-PCR/V

VOLUME: 2 ml

LOT: P071029-003

EXPIRATION: 6 August 2011

VIRUS: Porcine Circovirus (PCV2)



DESCRIPTION: Anti-PCV2 polyclonal antiserum. Liquid. Porcine origin. Made against a field isolate of type II; does not stain type I.

INTENDED USE: Useful for IFA. Has not been successfully applied to immunohistochemistry on formalin-fixed, paraffin-embedded tissues.

QUALITY CONTROL METHOD: Indirect FA using VMRD PCV2 12-well Substrate Slides (catalog no. SLD-IFA-PCV2); Anti-Pig IgG FITC Conjugate (catalog no. 085-1 or 085-10).

Specific Reaction: 4+ positive fluorescence at 1:1600 with no background. Endpoint titer is at 1:25,600.

Other Reactions or Comments: NA

STORAGE: This antiserum is provided in liquid form and should be stored at 4-8°C. **DO NOT FREEZE!** If antiserum becomes cloudy, it should be discarded. This antiserum contains 10 ppm ProClin 300 as a preservative.

FOR *IN VITRO* LABORATORY USE ONLY.

WARRANTY: VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.