


CERTIFICATE OF ANALYSIS

2E10-H2

Monoclonal Antibody

<b>Catalog No. / Cell Line:</b>	2E10-H2
<b>Lot:</b>	P080821-004
<b>Isotype:</b>	IgG <sub>1</sub>



**Specificity:**

Canine Adenovirus Type 1

**Known Applications:**

Virus neutralization.

**References:**

- Whetstone, C.A. Monoclonal antibodies to canine adenovirus 1 and 2 that are type-specific by virus neutralization and immunofluorescence. *Vet. Microbiol.* 16(1):1-8 (Jan. 1988).
- Whetstone, C.A. Should the criteria for species distinction in adenoviruses be reconsidered? Evidence from canine adenoviruses 1 and 2. *Intervirology* 23(2):116-120 (1985).

**Description and Handling:**

This monoclonal antibody is produced as mouse ascites fluid, clarified by centrifugation, and filtered through a 0.2  $\mu$ M filter. The concentration is 1.0 mg/ml in phosphate-buffered saline, preserved with sodium azide.

**Storage:**

When the vial is stored at 2-8°C, it should be stable for one year.

**FOR IN VITRO LABORATORY USE ONLY.**